

USER-MANUAL

TLC40 / TLC80

TLC80DP

TLC90



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 1/33



est. 1878

1	SAFETY AND WARNINGS.....	3
2	WARRANTY.....	3
3	PRECAUTIONS AND HAZARDS	4
4	INSTALLATION.....	4
4.1	UNPACKING	4
4.2	PLACEMENT AND VENTILATION.....	5
4.3	INITIAL USE	5
4.4	BATH LIQUIDS	6
4.5	FLUID LEVEL	7
4.6	MAXIMUM FLUID LEVEL.....	8
5	CONNECTING	9
5.1	CLEANING	9
6	INTRODUCTION TO THE TLC SERIES	10
6.1	GENERAL	10
6.2	PUMP	10
6.3	TEMPERATURE CONTROL	11
7	OPERATION	12
7.1	SWITCHING ON.....	12
7.2	CAUTION WITH POWERING ON/OFF	12
7.3	CONTROL PANEL.....	13
7.4	SELECTING PARAMETERS	13
7.5	QUICK START.....	14
7.6	SELECTING CAPACITY (TLC80DP ONLY)	15
7.7	TUNING THE BATH	15
7.8	TUNE PID VALUES.....	16
7.9	DESELECT TUNING	16
7.10	MANUAL TUNING	16
7.11	DRAINING BATH FLUID	17
7.11.1	Using the drain tap	18
7.11.2	Fast method: drain by hand.....	18
7.12	MAINTENANCE	19
8	TROUBLE SHOOTING	19
8.1	GENERAL	19
8.2	ERROR MESSAGES ON DISPLAY.....	21
8.3	FAULTY TEMPERATURE READINGS	21
9	SPECIFICATION.....	21
9.1	TECHNICAL SPECIFICATIONS OVERVIEW	21
9.2	TLC40 PERFORMANCE	23
9.2.1	Stability TLC40	24
9.3	COOL DOWN TLC40.....	24
9.4	TLC80 - TLC90 PERFORMANCE	25
9.4.1	Stability TLC80 - TLC90.....	25
9.4.2	Cool down TLC80 - TLC90	25
9.4.3	Heat removal TLC80 - TLC90.....	25
9.5	TLC80DP PERFORMANCE	26
9.5.1	Stability TLC80DP.....	26
9.5.2	Cooldown TLC80DP	26
9.5.3	Heat removal capacity TLC80DP	26
10	SPARE PARTS	27
10.1	TLC80 AND TLC90	27
10.2	TLC80 DP.....	27
10.3	TLC40	28
11	PRINCIPLE OF OPERATION	29
12	EC DECLARATION OF CONFORMITY.....	31
13	DISCLAIMER.....	33

ISO 9001 : 2015
NL/PRO 238239125van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 2/33

1 SAFETY AND WARNINGS

Make sure before installing or operating the equipment to read and understand all instructions and safety precautions listed in this manual. If there are any questions concerning the operation of the equipment or about the information given in this manual please contact your local dealer or our sales department first.

Performance of installation, operation, or maintenance other than those described in this manual may result in a hazardous situation and may void the manufacturer's warranty.

Never operate equipment that is not correctly installed. Unqualified personnel must not operate the equipment. Avoid damage to the equipment, or its accessories, caused by incorrect operation.

Important:

- When performing service, maintenance or moving the apparatus, always disconnect the line cord of the apparatus,
- Proper skilled and trained personnel are only allowed to operate this equipment,
- Take notice of warning labels and never remove them,
- Refer service and repairs to qualified technician,
- If a problem persists, call your supplier or Tamson Instruments B.V..

2 WARRANTY

Tamson Instruments B.V. warrants that all their manufactured equipment is free from defects in material and workmanship, preventing the machine from normal operation. Tamson Instruments B.V. does not warrant that the equipment is fit for any other use than stated in this manual. The manufacturer can only be held responsible for the security, reliability and performance of the equipment, when operated in accordance with the operating instructions, extensions, adjustments, changes and/or if repair is performed by Tamson Instruments B.V. or authorized persons only. This warranty is limited to one year from the date of invoicing. All equipment and materials are subject to standard production tolerances and variations.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 3/33

3 PRECAUTIONS AND HAZARDS

Before attempting to operate the bath read all parts of this manual carefully to insure smooth operation and avoid damage to the equipment or its accessories.

If a malfunction occurs, consult section "TROUBLE SHOOTING ", page 19

If problem persists, call your supplier or Tamson Instruments B.V. Never operate the equipment if not correctly installed. The equipment must be operated only by qualified personnel. Avoid damage to the equipment or its accessories through incorrect operation.



**READ
CAREFULLY**


4 INSTALLATION

Tamson Instruments B.V. is not responsible for any consequential damage or harm caused by using this bath. Repairs on the electrical system of the bath may only be carried out by well trained and authorized persons.

4.1 Unpacking

To avoid damage during transport all Tamson baths are carefully packed for shipment. Check the packaging for external damage and make a note on the shipping documents if any damage is found. Always retain the cartons and packing material until the bath has been tested and found in good condition. Transport companies generally will not honor a claim for damages if the respective box is not available for examination. The shipment contains at least the bath/circulator. The consignment may contain other parts, individually packed in small boxes. Please see packing list for details concerning total contents of consignment.

Before filling the bath remove any remaining packing material from its interior. The interior of the bath can be accessed by taking off the lid on the top of the bath.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370


NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 4/33


4.2 Placement and ventilation

The bath has to be placed in a well ventilated area.



Do not place the bath in a dusty environment. Dust will block the condenser unit inside the apparatus which will lead to severe mechanical damage.


Air circulation has to be enabled by 10 cm of free space at both sides and 15 to 20 cm at the rear of the bath.



4.3 Initial use

Due to transportation allow the bath a 24 hour period for stabilization. This applies also when the apparatus has been tilted or fallen.

Do not switch on power, because lubrication-oil inside the compressor system has to run into the capillary. It will take several hours before the oil has flown back into the compressor unit. When the apparatus is immediately turned on after it has been tilted, severe damage may occur to the compressor unit due to insufficient lubrication.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046


Chamber of commerce 27 16 95 41

Page 5/33

4.4 Bath liquids

The bath must be filled with a liquid suitable for the minimum operating temperature.

It is very important to select a liquid with a viscosity of less than 20 cSt at the operating temperature and a flash point which is well above the operating temperature.




The best bath liquid to be used is methanol, enabling circulation at very low temperatures.

Methanol is extremely flammable and can cause fire hazard. Take all necessary precautions to reduce fire hazard. When using methanol all vapors must be removed by using appropriate air ventilation.

Methanol is toxic and can cause health risks. Use appropriate ventilation and other precautions to prevent inhaling toxic vapors. If ventilation is insufficient the risk of explosion hazards can occur!

The supplier of the bath liquid (methanol) will be able to hand over all chemical details and safety precautions related to the use of methanol. These precautions must be followed when operating the bath.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 6/33

4.5 Fluid level

Do not operate the bath with low fluid level.

When the fluid level is too low, bath fluid will vaporize leading to toxic and flammable fumes.

Flammable fumes can lead to fire



Flammable fumes can be ignited bij the not submersed part of the heating element.

The level indicator(optional) will start to blink (blue light) when the fluid level is too low.

When the bath has been installed it must be filled with an appropriate liquid. When working with water the bath should be filled to 1 cm below the lid. For oil the bath should be filled to not more than 5 cm below the lid. Depending on the operating temperature the liquid level in the bath should be observed and excessive fluid should be removed.

The liquid level should be maintained between 1 and 3 cm below the lid during normal operation.

The heating element will be damaged when not fully submerged in the bath fluid. A lower level than 5 cm below the lid may damage the heaters. A high bath level can cause overflow and will might also damage the bath insulation.



Low fluid level

- Can cause fire when heater get's partially exposed
- Will damage the heater



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.



E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 7/33

4.6 Maximum fluid level

If the fluid level is too high, it will leave the bath via the overflow outlet (10mm outside diameter pipe). Prevent fluid from the overflow outlet entering the backside of the apparatus. For this reason the overflow outlet must be connected to a waste container.

When the bath is working at low temperatures, tubing and waste container must be chemical resistant and able to withstand high temperatures (<-80°C / -112°F).



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046


Chamber of commerce 27 16 95 41

Page 8/33

5 CONNECTING

Before plugging TLC bath into mains socket, make sure the voltage of the bath corresponds to the local voltage and frequency.




Use a mains supply that is well earthed, clean of interference and suitable for the acquired electrical load of the bath.



5.1 Cleaning


Regularly check the apparatus and condenser unit for dust. Follow steps 1 to 6 to remove the dust.

Unplug from mains supply



Remove dust with a vacuum cleaner

Do not use compressed air to clean



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370


NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 9/33

6 INTRODUCTION TO THE TLC SERIES

The TAMSON model TLC circulators are designed to perform accurate temperature control required for general laboratory. The Tamson circulator is optimized for temperature control of applications requiring a high degree of stability over a broad temperature range.



6.1 GENERAL

The TLC apparatus consists of a combination of a cooling system and a microprocessor controlled heating element. This design ensures a high degree of accuracy and reproducibility of temperature controls.


The TAMSON baths are constructed throughout from corrosion-resistant – stainless steel and Teflon – materials. The bath is effectively insulated against heat loss by a layer of Armaflex® rubber between the inner tank and outer casing.

6.2 Pump

A circulation pump is built-in to guarantee an uniform temperature distribution within the bath and providing the possibility to circulate through a closed external system.

The TLC40, TLC80, 80DP and TLC90 can be provided with two pumps.

Standard pump: a circulation pump is provided with a circulation pump with a capacity of 10 liters / min and a max pressure of 3 m head of water [3mtr H₂O].



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.


E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 10/33

A performance pump: can be ordered under item 24T0392. This pump is not interchangeable and has to be specified when ordering as an optional item.


The standard pump offers a motor protection (small thermal fuse) which is placed on the back panel of the motor compartment. This protects the motor from excessive loads.

When pumping viscous liquids the thermal fuse can be activated. The motor fuse can be reset by pressing the knob.



The inlet and outlet of the pump are fitted with 10 mm hose connections and 3/8" threaded fittings for metal tubing. The outlet (most near to the backside of the apparatus) is provided with a small disc, stopping the circulation action.


Before the fluid can begin circulating through external systems this disc must be removed by unscrewing the fitting and hose connection removing the disc.



Optional available is a stronger pump with adjustable flow. Part number is 24T0392. Flows of up to 15 liters per minute can be realized. This pump is for circulating purposes where a constant external flow with higher pressure up to 0.8 bar is required. Please check our specification sheet for all pump performance details.

6.3 Temperature control

Temperature control and setting the bath temperature is regulated using a PT-100 temperature probe Class A connected to a microprocessor module. The advanced electronic control system continually computes the energy input required for optimal temperature accuracy and stability. Temperature read-out is on a 4 digit, 7 segments LED-display with 0.1°C reading. The actual set point is computed within the controller with an accuracy of 0.01°C.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 11/33

7 OPERATION


7.1 Switching on

If the bath has been properly filled with fluid it can be switched with the mains switch located on the front panel. Choose a working temperature (set point) with **▲** and **▼** keys. For an optimal performance the bath is best tuned at the set point temperature.

7.2 Caution with powering on/off

Be careful and do not toggle with the on/off switch. To start the cooling compressor high currents are needed which will heat the compressor motor internally. Switching off and on the compressor several times will lead to mechanical damage.

When switching off the apparatus, wait for 10 minutes before switching the system back on again.



1x



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370


NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 12/33

7.3 Control panel

- 1 LED display
- 2 Heater indicators (Green)
- 3 Level indicators (Blue)
- 4 Over temperature indicator (Red)
- 5 Up and down buttons
- 6 List – key
- 7 Page - key



7.4 Selecting parameters

The front panel layout shows the following 4 operating keys:



Press shortly

- | | |
|------|---|
| SP.1 | - Temperature set point in °C, |
| tU.n | - Tuning the bath ("Atune"), |
| Pb | - Changing the tuning (PID) parameters "P" value |
| ti | - Changing the tuning (PID) parameters "I" value |
| td | - Changing the tuning (PID) parameters "D" value |
| OFSt | - Change temperature readout offset ($\pm 5^{\circ}\text{C}$) |

Sequence

Each time

Up and Down keys allow changing the listed value. All changed values, like set point and PID parameters, will be kept in memory when pressing the page key to leave the menu. After switching off the power supply, changed values are kept in memory.

- SP.1 - Set point adjust. Press page key to confirm setting when altered.
- tU.n - Auto tuning PID parameters, set on or off.
Off is default value here.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 13/33

On, when selecting parameter 2nd time, On.A is played. This indicates adaptive tune. Bath is adapting PID setting in background operation. When in On.A mode PID values can not be set manually anymore.

- pb^d** - Changing the tuning (PID) parameters. 'Proportional band value. Default value is 2, range 1 .. 100. (Proportional band equals 1/proportional value). Press page key to confirm setting when altered.
- ti^d** - Changing the tuning (PID) parameters 'Integrator' value. Default value is 200. Range 1 .. 1200. Default value is 200. Press page key to confirm setting when altered.
- td^d** - Changing the tuning (PID) parameters 'Differentiator' value. Default value is 10. Range 1 .. 600. Default value is 40. Press page key to confirm setting when altered.
- OFST** - Set offset. The offset value is added to the displayed bath temperature. I.e. when the displayed bath temperature is 20.0°C, an offset of 0.2 displays a value of 20.2°C. An offset of -0.2°C will display 19.8°C. Range is -5°C .. +5°C. Press page key to confirm setting when altered.

note:

- ☞ When tune is set to on, PID values can not be set. Change tune to off en alter PID values.
- ☞ Press page key to confirm setting when altered.

7.5 Quick start

To quickly start operating the bath:

- Fill the bath with sufficient fluid* to flood the cooling coil, However do not fill the bath fully to the maximum.
- Place the power plug,
- Switch the bath on using the mains switch,
- Choose a working temperature (set point):

Press  once. Display will indicate "SP.1". Alter set point temperature by using up  and down  key to select the desired bath temperature.

Press  to display bath temperature again (or wait 10 seconds).



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046


Chamber of commerce 27 16 95 41

Page 14/33

7.6 Selecting capacity (TLC80DP only)

The TLC80DP system is constructed around 2 x 2 compressors. The system contains a first and second stage and each stage has two compressors. In some cases the full heat removal power is not required(idle or overnight use). In these circumstances each second compressor in both stages can be switched off. The purpose for this is to conserve energy and when heating up the bath prevent high pressure in the cooling circuit.

When the bath is working at low temperatures and a new much higher temperature is entered as SP, always switch back the TLC80DP to half capacity.



7.7 Tuning the bath

The temperature control of the bath is based on a digital PID system. When using different fluids in the bath each with their own heat capacity, the use of external cooling and external connected processes (circulation), or working at different set point temperatures requires new settings of the PID parameters. These parameters have to be optimized after changes to the system when optimal and accurate temperature control of the bath liquid is required.

Tuning of the bath results in:

- Stable temperature control of the bath,
- No over- or undershoot of the temperature set point,
- Quick response to deviations from the set point caused by external disturbances.



Tuning can be done automatically or manually. The parameters mentioned in Table 1 will influence the control of the bath.

Parameter	Description	Display	Default value
Proportional band	The bandwidth in display-units over which the output power is proportional between minimum and maximum	Pb	2
Integral time	Determines the time taken by the controller to remove steady state error signals	Ti	40 .. 200
Derivative time	Determines the time taken by the controller to react on error signals	Td	0



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 15/33

7.8 Tune PID values

Choose desired set point,
Wait until the bath has reached the set point temperature,
Highest temperature accuracy can be achieved when
tuning the bath. Select tU.n by pressing twice.

Press to choose "ON",

Press 5 times to return to displayed PV
(PID values can no longer be adjusted unless the "tU.n"
value is set to off)

7.9 Deselect tuning

- Pressing twice,
- Display shows "tU.n",
- Press to choose "Off",
- Press 5 times to return to displayed PV.

7.10 Manual tuning

The parameters for the PID control can also be changed manually. The method described below provides fast finding of the PID settings and is referred to as the method of "Ziegler Nichols"

This procedure is only to be followed when automatic tuning is not functioning.



This method provides fast manual findings of the PID values. Start the bath at its required running temperature:

1. Set the integral time "Ti" and the derivative Time "Td" to off
2. Check if the Lcb and Hcb are set to auto
3. Ignore the fact that the temperature may not settle precisely at the set point
4. If the temperature is stable, reduce the proportional band Pb so that the temperature just starts to oscillate. If the temperature is already oscillating, increase the proportional band until it begins oscillating. Allow enough time between each adjustment for the loop to stabilize. Make a note of the proportional band value "B" and the period of oscillation "T". Set the Pb, Ti and Td parameter values according to the calculations given in the next table.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046


Chamber of commerce 27 16 95 41

Page 16/33

Type of control	Proportional band "Pb"	Integral time "Ti"	Derivative time "Td"
Proportional only	2*B	Off	Off
P + I control	2,2*B	0,8*T	Off
P + I + D control	1,7*B	0,5*T	0,12*T


7.11 Draining bath fluid

Before removing flammable bath liquids take the appropriate fire hazard precautions against these liquids.



Do not remove cold bath liquid. Cold liquid can cause severe burning when spilled. Preheat the bath liquid back to ambient temperature before removing.

When removing bath fluid do not inhale toxic vapor. Always use appropriate ventilation.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com


VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 17/33

7.11.1 Using the drain tap

The TLC baths can be drained via the drain tap located at the backside of the apparatus.


For safety reasons the tap can only be opened by using a screwdriver. The thread inside the tap is 5/8".

**USE SCREWDRIVER
TO OPEN TAP
5/8" BSP inside thread**



7.11.2 Fast method: drain by hand

The bath fluid can easily be removed by using a length of hose. Length of approximately 1.5 mtrs and 10mm inner diameter.




A.


Handle as follows:

Fill the hose with fluid ("A"),

- Close one side of the hose and immerse the other end in the bath ("B"),
- Lower the closed end into a waste bucket and let the bath fluid flow ("C"),
- Remove the last pieces with a sponge or tissue,
- Refill the bath approximately 1 cm below the lid.



B.



C.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 18/33

7.12 Maintenance

Keep the apparatus free from dust. Regularly check the cooling openings and remove dust with vacuum cleaner. If necessary remove cover and clean internally. Use appropriate protection when cleaning, dust can be very unhealthy.

8 TROUBLE SHOOTING

8.1 General

The motor is not running and the lamps are not on.	Check the main voltage and all electrical connections, including switches. Check fuse in wall socket of mains supply.
Motor is not running, lamps and electronics operate (the motor turns freely by hand).	Motor overloaded. The motor protection may trip because of mechanical overloading of the pump. Reset fuse by pushing button on rear side of blue motor housing.
	The motor capacitor has a defect. Replace motor capacitor inside the system.
Bath is not heating, Red error led on front panel blinks on and off	Malfunctioning wiring. Check wiring. Over temperature protection has been activated. Bath temperature reached > 30'C. Switch off the unit, Let the bath fluid cool down Restart.
Compressor is making "clicking" noise	Bath has been switched on and off too quickly. Wait approximately 10 minutes before switching on again. The clicking noise is a temperature fuse inside the compressor preventing it from overheating.
	Bath liquid temperature is rising too fast. When the bath is working at low temperatures and a new much higher temperature is entered as SP, always switch back the TLC80DP to half capacity. This will speed up the heating of the bath liquid but more important it will prevent excessive high pressures in the cooling circuit.
Compressor will not start	When the machine has been switched off the pressure in the cooling circuit can rise to a high level. For this reason the apparatus has to settle the pressure of the cooling liquid for several minutes (5 ... 10) before it can be restarted.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 19/33

When this problem occurs more often it will indicate that dust is blocking the condenser circuit. Clean the inside of the TLC80 and remove all dust from the condenser unit and fans with a vacuum cleaner. Do not use pressurised air to remove dust. The flying dust particles can damage the fan bearings, are very unhealthy when breathed in and can lead to severe eye problems.

Can not set PID values

Check if ATUN has been set to on. PID values can no longer be adjusted unless the ATUN value is set to off. See 7.7 page 15.

Temperature not stable

Tune the TLC80 at the set point temperature. Place top lid on the bath.

Temperature does not drop

Bath is not properly insulated. Place top lid on the bath.

Insulate circulation circuit.

Circulation load exceeds heat removal capacity.

Use appropriate bath fluid preferably methanol. If the bath contains a high percentage of water, this will freeze and prevent the fluid from cooling, as can be seen from the picture on the left.



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 20/33

9.2 Error messages on display

Code	What does it mean	What action to take
xx	Error on configuration parameter	The 2 digits show the wrong parameter.
101	Error on RTD	Replace PT100
102	Error on voltage input calibration	Defective controller, replace
130	Error on operative mode parameter	Contact Tamson or local dealer
150	Error during data storing	Fault resets itself
Display "dead"	Probably over-temperature protection is active	Turn bath off and on
Red LED	Blinking	Over temperature protection has been activated. Bath temperature reached > 30°C. Switch off the unit, Let the bath fluid cool down Restart.

8.3 Faulty temperature readings

The temperature read out on the display does not correspond to the temperature measured.

PT100 is defective,
Offset value is set to a wrong value. See
Selecting parameters page 13, parameter
"OFST".

9 SPECIFICATION

9.1 Technical specifications overview

Item	Unit	TLC15	TLC30	TLC40	TLC80	TLC80DP	TLC90
Temperature							
Range	[°C]	+60..-15	+60..-30	+20..-40	+20..-80	+20..-80	+20..-90
Setting ±	[°C]	0,1	0,1	0,1	0,1	0,1	0,1
Stability* ±	[°C]	0,05	0,05	0,05	>0,05	>0,05	>0,05
Uniformity* ±	[°C]	0,05	0,05	0,05	>0,05	>0,05	>0,05
Specified at		-5°C)	-5°C)	-5°C)	-5°C)	-5°C)	-5°C)
Heating							
Rate heating	[°C/min]			2,3	2,3	2,3	
Heater	[W]			1100	1500	1500	1500
Pump							
Pressure	[mbar]			300 max	300 max	300 max	300 max
Capacity	[L/min]			10	10	10	10
Pump high pressure							
Pressure				>800	>800	>800	>800
Capacity				15	15	15	15
Dimensions							
Bath volume	[L]			14..15	14..15	14..15	14..15



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 21/33

Bath Opening	[mm]			240 x 170	240 x 170	240 x 170	240 x 170
Opening effective				240 x 160	240 x 160	240 x 160	240 x 160
Depth	[mm]			150	150	150	150
Length	[mm]			810	810	789	810
Width	[mm]			460	460	496	460
Height	[mm]			770 (710 no wheels)	770 (710 no wheels)	1175	770 (710 no wheels)
Weight	[kg]			80	80	140	80
Power (max)**	[kW]			2.4 max	2.8 max.	3.6 max	2.8 max.
Voltage	[Volt]			230 or 115	230 or 115	230	230 or 115

* Absolute min/max value measured over 1hrs in methanol

** Depends on bath temperature and cooling or heating cycle



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:


NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 22/33

9.2 TLC40 performance



ISO 9001 : 2015
NL/PRO 238239125


van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.


E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 23/33

9.2.1 Stability TLC40



9.3 Cool down TLC40



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117


Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com


VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 24/33

9.4 TLC80 - TLC90 performance

9.4.1 Stability TLC80 - TLC90



9.4.2 Cool down TLC80 - TLC90



9.4.3 Heat removal TLC80 - TLC90

Temperature [°C]	TLC80	TLC90	
20*	500	500	[watt]
0*	450	450	[watt]
-20*	350	350	[watt]
-50*	250	250	[watt]
-60*	225	225	[watt]
-70	175	175	[watt]
-80	125	125	[watt]
-90	-	20	[watt]

* During cool down. Do not operate bath continuously at this temperature



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370


NL95 RABO 0160100046

Chamber of commerce 27 16 95 41


Page 25/33

9.5 TLC80DP performance

9.5.1 Stability TLC80DP



9.5.2 Cooldown TLC80DP



9.5.3 Heat removal capacity TLC80DP

Temperature [°C]	TLC80DP	[watt]
20*	550	
0*	500	
-20*	450	
-50*	400	
-60	400	
-70	325	
-80	200	

* During cool down. Do not operate bath continuously at this temperature



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 26/33

10 SPARE PARTS

10.1 TLC80 and TLC90

230 Volt 50HZ	230V 60HZ	115 Volt 60Hz	Description
25T1290	25T1300		Motor for pump
24T3300	n. a.		Capacitor 7uF
n. a.	24T3330		Capacitor 25uF
25T0253	25T0254		Heater 1100 Watts
24T8081	n. a.		Motor fuse 0.3 Amp.
n. a.	24T8080		Motor fuse 0, 6 Amp.
24T8581			Mechanical safety thermostat
04T2070			Bearing cover with bearing
04T2135			Bearing PTFE
22T3520			Fill valve
24T8545			Mains switch
28T4015			Front foil
25T2310			PT-100 sensor
06T0475			Mains board with filter and relay
28T3015			Controller

10.2 TLC80 DP

230 Volt 50HZ	Description
25T1295	Motor for pump
25T1343	Capacitor 2uF
25T0251	Heater 1100 Watts
24T8081	Motor fuse 0.3 Amp.
24T8581	Mechanical safety thermostat
04T2070	Bearing cover with bearing
04T2135	Bearing PTFE
24T8545	Mains switch
28T4015	Front foil
25T2310	PT-100 sensor
06T0475	Mains board with filter and relay
28T3015	Controller



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 27/33

10.3 TLC40

230 Volt 50HZ	230V 60HZ	115 Volt 60Hz	Description
25T1290	25T1300		Motor for pump
24T3300	n.a.		Capacitor 7uF
n.a.	24T3330		Capacitor 25uF
25T0251	25T0252		Heater 1100 Watts
24T8081	n.a.		Motor fuse 0.3 Amp.
n.a.	24T8080		Motor fuse 0,6 Amp.
06T0465	06T0470		PCB with relay and mains
24T8581			Mechanical safety thermostat
04T2070			Bearing cover with bearing
04T2135			Bearing PTFE
22T3520			Filling valve
24T8545			Mains switch
28T4015			Front foil
25T2310			PT-100 sensor
28T3014			Controller



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370


NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 28/33

11 PRINCIPLE OF OPERATION

TLC80, TLC80DP and TLC90: principle of operation			
1	Compressor(s) 2nd stage	10	Compressor(s) 1st stage
2	Safety overpressure switch 2nd stage	11	Safety overpressure switch 1st stage
3	Oil separation (DP only)	12	Heating top plate
4	Inter cooler (DP only)	13	Fan condenser
5	High pressure	14	Condenser
6	Heat exchanger	15	High pressure side
7	Throttle second stage	16	Dryer
8	Low pressure side 2nd stage	17	Throttle 1st stage
9	Low pressure side 1st stage		



ISO 9001 : 2015
NL/PRO 238239125


van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 29/33

TLC40			
1	Compressor	5	Filter / dryer
2	Schräder - High pressure side	6	Thermostatic Expansion Valve
3	Fan condenser	7	Evaporation (Bath)
4	Condensor	8	Schräder - Low pressure side (suction)



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 30/33

12 EC DECLARATION OF CONFORMITY

Following equipment is in compliance with EMC Directive 2014/30/EU:



Product: Thermostatic bath and circulator
 Model: TLC40
 TLC80
 TLC80 DP
 TLC90
 Serial code: Effective from 02Txxx
 Manufacturer: Tamson Instruments bv
 van 't Hoffstraat 12
 2665 JL Bleiswijk
 The Netherlands

The products are in conformity with the following specifications:

Item	Reference	Description	Test result
a	RoHS Directive	2011/65EU	p
b	EN61010-2-010	Safety requirements for electrical equipment for measurement, control, and laboratory use. Particular requirements for laboratory equipment for the heating of material	
c	Machine Directive 2006/42/EC	Machinery Directive, of the European Parliament and of the Council of 17 May 2006/42/EC 2nd Edition June 2010	p
d	EN 60204	Machinery Directive and Safety requirements	p, p ⁱ
e	EN60950-1	Low Voltage Directive	p
f	EN61000-3-2:2014	Harmonics	p
g	EN61000-3-3	Flicker	p ³
h	EN61000-4-2 +A1+A2	ESD	p
j	EN61000-4-3 +A1+A2	Radiated immunity	p (anechoic room)
k	EN61000-4-4	Electrical Fast Transients	Minimum requirements pass
l	EN61000-4-5+A1	Surges	Minimum requirements pass
m	EN61000-4-6+A1	Conducted immunity	p
n	EN61000-4-11 +A1	Voltage dips and Voltage variations	p
o	EN55016-2-1	Conducted emission	p
p	EN55016-2-3	Radiated emission	p (anechoic room)
q	Pr EN 378	Refrigerating systems and heat pumps - Safety and environmental requirements	
r	EN 13445-5	PED Inspection and Testing	Maximum working pressure level of 20 Bar is confirmed. On each apparatus following pressure and leak tests have been carried out with positive result - Low pressure side 10 Bar - High pressure side 25 Bar

p = Pass

pⁱ = Individually tested

p³ = Pass, condition of operating during Pst measurement: Operational with heating element 1400W.
 P_{st} and P_{lt} are not evaluated in accordance with A.5 of Annex A of EN 61000-3-3(1995) + A1(2001).



ISO 9001 : 2015
 NL/PRO 238239125

van 't Hoffstraat 12
 2665 JL Bleiswijk, The Netherlands
 T. 31 (0) 10 522 43 73
 tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
 Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 31/33

not applicable were:

Conducted discontinuous emissions (Clicks)

Radiated emission (OATS)

Magnetic field immunity

The equipment conforms with all the specifications and norms in this regard.

The equipment conforms without any further notice.

Entity responsible for marking this declaration :

Manufacturer, Tamson Instruments bv, van 't Hoffstraat 12, Bleiswijk The Netherlands,

Name	:	R.C. van Hall
Function	:	Director
Date	:	January, 2018
Version	:	1.04



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01

Bank account no.:

NL28 INGB 0007 350 370

NL95 RABO 0160100046

Chamber of commerce 27 16 95 41

Page 32/33

13 DISCLAIMER

Information given herein is offered in good faith as accurate, but without guarantee. Conditions of use and suitability of the product for particular uses are beyond our control; all risks of use of the product are therefore assumed by the user.

We expressly disclaim all warranties of every kind and nature, including warranties of merchantability and fitness for a particular purpose in respect to the use or suitability of the product.

Nothing is intended as a recommendation for uses which infringe valid patents or as extending license under valid patents. Appropriate warnings and safe handling procedures should be provided to handlers and users. Alteration of this document is strictly prohibited. Except to the extent required by law, republication or retransmission of this document, in whole or in part, is not permitted. Tamson Instruments bv assume no responsibility for accuracy of information. Tamson Instruments neither represent nor warrant that the format, content or product formulas contained in this document comply with the laws of any other country except the Netherlands.

Tamson Instruments bv, All rights reserved



ISO 9001 : 2015
NL/PRO 238239125

van 't Hoffstraat 12
2665 JL Bleiswijk, The Netherlands
T. 31 (0) 10 522 43 73
tlc40-80-80dp-90man Rev. 3.01 UK 1117

Tamson Instruments B.V.

E-mail: info@tamson.com
Website: www.tamson.com

VAT: NL 80 66 34 984 B01
Bank account no.:
NL28 INGB 0007 350 370
NL95 RABO 0160100046
Chamber of commerce 27 16 95 41
Page 33/33