

ACCESSORIES

Manifold

REF 00T0635

APPLICATION

To connect two or more devices to one cold flow of a Tamson circulator.

MAIN CHARACTERISTICS

The manifold unit can be mounted on the wall and uses manual valves connecting the analysers to the circulating cold flow. The manifold can be used to connect 2 devices to a single low temperature circulator. The manifold can be set in series offering 4 devices to be connected to one cooling circulator. Its small dimensions, 270 x 126 mm (WxH), makes it easy to mount and install.

The manifold is well insulated. Tamson also provides low temperature insulated hoses to safely and easy connect the cold property apparatus. These hoses consist of a flexible stainless steel inner hose and are insulated with black armaflex covered with a protecting nylon sleeve. All connections are standard 3/8" BSP. For interconnection to other threads, stainless steel adaptor pieces are available.

